

# golden **COMPASS™** **HD**

## **Manual de usuario**



***Compass HD™***  
**Silla mecánica de tracción central**  
**GP620 SS/GP620 CC**



401 Bridge Street, Old Forge PA 18518  
Tel: (800) 624-6374 Fax: (800) 628-5165  
[www.goldentech.com](http://www.goldentech.com)



Gracias por comprar su silla mecánica *Compass HD GP620*. Nos sentimos honrados de que haya elegido a Golden Technologies para sus necesidades de movilidad.

*Compass HD* Modelo GP620

- Incluido en el paquete:
- 1 Manual de usuario
  - Unidad base
  - 1 Asiento
  - 2 reposabrazos
  - 1 Banderilla de seguridad
  - 1 Mando de control
  - 1 Cinturón de seguridad
  - 1 Cargador de batería

**PARA SU REFERENCIA**

Por favor llene los datos de su *Compass HD™*. Esta información será útil en el caso que alguna vez necesite contactar a Golden Technologies, Inc. en relación a su silla mecánica.

**Su *Compass HD™***

Modelo \_\_\_\_\_ Número de serial \_\_\_\_\_

Fecha de compra \_\_\_\_\_ Color \_\_\_\_\_

Opciones \_\_\_\_\_

**Su representante o proveedor de Golden Technologies, Inc.**

Nombre \_\_\_\_\_

Compañía \_\_\_\_\_

Dirección \_\_\_\_\_

**Por favor recuerde llenar y devolver la tarjeta de registro de su garantía**



**DESCARGO DE RESPONSABILIDAD**

Le pedimos que lea este manual completamente antes de operar su nuevo *Compass Sport™*. Golden Technologies, Inc. no es y no puede ser responsabilizado por ningún daño o lesión incurrida debido al uso inapropiado o inseguro del *Golden Compass Sport™*. Golden Technologies, Inc. específicamente renuncia a la responsabilidad de cualquier lesión corporal o daño material que podría ocurrir durante cualquier uso que no cumpla con las ordenanzas o leyes aplicables locales, estatales o federales.

**Uso previsto:** El uso previsto de la silla mecánica Compass HD es proveer movilidad para adultos limitados a una posición sentada y capaces de operar un simple control manual.

## CONTENIDO

INTRODUCCIÓN .....	4-5
SEGURIDAD.....	6-9
EMI/IRF.....	10
ESPECIFICACIONES.....	11
ARMADO.....	12-16
AJUSTES DE SEGURIDAD.....	17-22
OPERACIÓN/DIAGNÓSTICOS .....	23-29
CARGANDO LAS BATERÍAS .....	30-31
CUIDADO Y MANTENIMIENTO .....	32
DESMONTAJE.....	33
GARANTÍA .....	34-35
TARJETA DE REGISTRO .....	36





# SAFETY

Your *Compass HD™* is a battery operated personal mobility vehicle. Please exercise caution and consideration when you are operating it. Driving your *Compass HD™* carefully and thoughtfully will help ensure your personal safety and the safety of other



**MANDATORY! Read and follow the owner's manual before operating your power chair.**

**NOTE:** Before learning to operate your *Compass HD™*, have your Golden Technologies, Inc. representative determine if it is advisable for you to practice getting on and off your *Compass HD™* and operating it in the presence of an attendant.

## **BEFORE GETTING ON YOUR *Compass HD™***

- ◆ Check to be certain that the power is turned off. See “Operation” on pages 23-28.
- ◆ This will eliminate the possibility of accidentally activating the joystick and causing injury to yourself or to others.
- ◆ Check to be certain that your *Compass HD™* is not in freewheel mode. See the “Operation” section on pages 23-28.
- ◆ Flip up the armrests.
- ◆ If your chair is equipped with leg rests, move them to the sides.
- ◆ If your chair is equipped with a footrest, flip it to the “up” position.

## **GETTING ON YOUR *Compass HD™***

- ◆ While being careful not to entangle yourself in the leg rests or the footrest, seat yourself comfortably and securely on the seat.
- ◆ Adjust the leg rests or flip down the footrest.
- ◆ Flip down or replace the armrests.
- ◆ Fasten the lap belt, if your chair is so equipped.

## **GETTING OFF YOUR *Compass HD™***

- ◆ Make certain that the power is turned off.
- ◆ Unfasten the lap belt, if your chair is so equipped.
- ◆ Flip up the armrests.
- ◆ Flip up the footrest, or move the leg rests to the sides.
- ◆ While being careful not to entangle yourself in the leg rests or the footrest, carefully stand and step away from the chair.

## **MAXIMUM WEIGHT**

Your *Compass HD™* has been rated to a maximum payload (passenger and anything else being carried on the *Compass HD™*) of **450 pounds**. Exceeding the maximum weight rating will void the warranty.



**WARNING**

**Exceeding the maximum weight capacity will void your warranty and may result in injury to yourself and/or others.**

# SAFETY

## DRIVING ON AN INCLINE

- ◆ Drive with **caution** when attempting to negotiate any incline, even handicap access ramps.
- ◆ Try to keep your *Compass HD™* **moving** when climbing an incline. If you do come to a stop, restart and accelerate slowly and carefully.



### WARNING

- ◆ **Always** climb or descend a gradient by driving straight up or straight down the face of the slope.
- ◆ **Do not** traverse or drive across the face of a gradient.
- ◆ **Do not** attempt to negotiate an incline that is covered with snow, ice, cut or wet grass, leaves, or any other potentially hazardous material.
- ◆ **Do not** back down an incline.
- ◆ **Do not** try to descend or climb an incline that is greater than the recommended maximum incline of **6 degrees**.



### WARNING

**If, while you are driving down a slope, your chair starts to move faster than you feel is safe, gently release the joystick and allow your *Compass HD™* to come to a stop. When you feel that you again have control of your power chair, adjust the speed control to a lower setting then push the joystick forward and continue safely down the remainder of the slope.**

## MEDICATION

Always check with your physician to determine if any of the medications you are taking may affect your judgment and/or your ability to operate your *Compass HD™*.



### WARNING

**Do not connect or allow anyone except an authorized Golden Technologies, Inc. representative to connect any electrical or mechanical device to your *Compass HD™*. Unauthorized accessories will void the warranty and may cause injury.**

# SAFETY

## SAFETY RULES



### WARNING

- ◆ **Do not** attempt to use your *Compass HD™* on an escalator. Always use an elevator.
- ◆ **Do not** carry passengers on your power chair.
- ◆ **Do not** operate your *Compass HD™* if it is not functioning properly.
- ◆ **Use caution** when driving on soft or uneven surfaces such as grass, gravel, and on decks where there is no railing.
- ◆ **Never** drive on the roadway, except when you must cross the street.
- ◆ **Always** cross streets at intersections and use the most direct route, making sure that your path is clear and that you are visible to motor traffic.
- ◆ It is **not recommended** to drive your *Compass HD™* up or down a step or curb that is higher than 1 1/2 inches.
- ◆ **Never** back up or down a step or curb.
- ◆ **Never** operate your *Compass HD™* while you are under the influence of alcohol.
- ◆ **Do not** operate or store your power chair where it will be exposed to rain, snow, mist, and below freezing temperatures.
- ◆ **Do not** operate your power chair on slippery, icy or salted surfaces.
- ◆ **Never** sit on your power chair when it is in freewheel mode and on an incline or decline.
- ◆ **Do not** modify your power chair in any way that is not authorized by Golden Technologies.
- ◆ **Do not** disassemble the tire. If disassembly is required, have your authorized Golden Compass HD provider perform any necessary maintenance or repair.
- ◆ **Do not attempt to** inflate the tires of your Golden *Compass HD™*. Your Golden *Compass HD™* is equipped with foam filled flat free tires that do not require inflation.
- ◆ **Never** sit on your *Compass HD* when it is being transported.
- ◆ **Always** fasten down your *Compass HD* securely with an approved tie-down system while transporting your power chair.
- ◆ Use **caution** when operating your *LiteRider Envy* while wearing long baggie clothing or clothing with loose belts, which could become entangled in the wheels of the power chair.
- ◆ **Always** keep your feet firmly on the footplate surface when operating your *LiteRider Envy* to eliminate any potential of your feet becoming entangled in the wheels of the power chair.

# SAFETY

Please be sure to follow these important warnings for the Golden Compass HD power wheelchair:



**WARNING**

**NEVER TRANSFER ONTO OR OFF OF THIS POWER WHEELCHAIR USING THE SEAT BACKREST FOR SUPPORT DURING TRANSFER. THE SEAT BACK MAY FOLD DOWN AND MAY CAUSE YOU TO LOSE YOUR BALANCE AND COULD RESULT IN PERSONAL INJURY.**

Do not modify your power chair in any way that is not authorized by Golden Technologies, Inc.



**WARNING**

**DO NOT CONNECT OR ALLOW ANYONE EXCEPT AN AUTHORIZED GOLDEN TECHNOLOGIES REPRESENTATIVE TO CONNECT ANY ELECTRICAL OR MECHANICAL DEVICE TO YOUR ALANTE SPORT. UNAUTHORIZED ACCESSORIES WILL VOID THE WARRANTY AND MAY CAUSE INJURY.**

Please be sure to follow this important warning if you encounter a problem with your mobility device.



**WARNING**

**Should you encounter a problem with the operation of your mobility device such as a mechanical or electrical problem immediately stop using the product. Disconnect the batteries then contact the provider from whom you purchased your mobility device for evaluation of the problem and repair if needed.**

# EMI/RFI

The rapid development of electronics, especially in the area of communications, has saturated our environment with electromagnetic (radio) waves that are emitted by television transmitters, cellular phones, citizen's band radios (CBs), amateur radios (ham radios), wireless computer links, microwave transmitters, paging transmitters, etc. These electromagnetic (EM) waves are invisible and increase in strength the closer one gets to the source of transmission. When these energy waves act upon electrical devices and cause them to malfunction or to function in an erratic or uncontrolled manner, they are referred to as *Electromagnetic Interference* (EMI) or *Radio Frequency Interference* (RFI).

## **EMI/RFI AND YOUR *Compass HD*<sup>TM</sup>**

All electrically powered vehicles, including power chairs are susceptible to EMI/RFI. This interference could result in abnormal, unintended movement of a power wheelchair.



**WARNING**

**Unintended movement or brake release could cause an accident or injury.**

The FDA has determined that each make and model of power chair can resist EMI/RFI to a certain level. The higher the level of immunity, the greater the degree of protection from EMI/RFI measured in volts per meter (V/m). The FDA has also determined that current technology is capable of providing 20 V/m of immunity to EMI/RFI, which would provide useful protection against common sources of interference. This product has been tested and has passed an immunity level of 20 V/m.

## **EMI/RFI RECOMMENDATIONS**



**PROHIBITED!** Do not turn on or use hand-held personal electronic communication devices such as cellular phones, walkie-talkies, and CB radios while your power chair is turned on.

- ◆ Be aware of any nearby transmitters (radio, television, microwave, etc.) on your intended route and avoid operating your chair close to any of those transmitters.
- ◆ Turn off the power if your *Compass HD*<sup>TM</sup> is going to be in a stationary position for any length of time.
- ◆ Be aware that adding accessories or components or modifying your power chair may make it more susceptible to EMI/RFI.
- ◆ If unintended movement or brake release occurs, turn your power chair off as soon as it is safe to do so.
- ◆ Report all incidents of unintended movement or brake failure to your Golden representative or to Golden Technologies, Inc. at (800) 624-6374.



**WARNING**

**Turn off your power chair as soon as it is safely possible if unintended or uncontrollable motion occurs or if unintended brake release occurs.**

# SPECIFICATIONS

<b>Model</b>	<b>GP620 SS</b>	<b>GP620 CC</b>
Medicare Code	K0824	K0825
Maximum Incline	6 degrees	6 degrees
Weight Capacity	450 lbs.	450lbs.
Maximum Speed	4.5 mph	4.5 mph
Operating Range <sup>1</sup>	21 miles	21 miles
Ground Clearance	3.5"	3.5"
Turning Radius	20 "	20 "
Type Batteries	2-NF22	2-NF22
Colors	Red, Blue	Red, Blue
Charger	Off Board	Off Board
Controller	60 Amp	60 Amp
<b>Dimensions:</b>		
Length	40"	40"
Width (outside of tires)	25"	25"
Height (ground to top of back of seat)	38.25" – 42.75"	38.25" – 42.75"
Height (ground to top of back of headrest)	47.25" – 53.75"	47.25" – 53.75"
Ground to Top of Seat	19.75" – 22.25"	21.25" – 24.25"
Number of Seat Height Positions	5	5
Size of Increments	0.5"	0.5"
<b>Weight of Unit: (assembled)</b>		
Base	101 lbs	101 lbs
Seat Arms	14 lbs	14 lbs
Seat (18 X 18 Solid Pan Seat)	34 lbs	
Seat (18 X 18 Captain's Chair Seat)		37 lbs
Batteries	NF22 - 37 lbs	NF22 - 37 lbs
<b>Tires:</b>		
	Flat Free	Flat Free
Center Drive Wheels (2)	10"	10"
Casters (4)	6"	6"
<b>Standard Seat:</b>		
	Captain's Van/Pan	Captain's Highback
Back Height-no headrest	23"	21"
Back Height with headrest	30.5" – 34.5"	28.5" – 32.5"
Width x Depth	20" x 20"	20" x 20", 22" x 20", 24" x 20"
Color	Black Vinyl	Black Vinyl
Electric Rise Seat (Option)	N/A	N/A
<b>Footrest:</b>		
Height adjustable	3" - 7"	3" - 7"
Angle adjustable	0 – 35 degrees	0 – 35 degrees
<b>Warranty:</b>		
Frame	5 Years	5 Years
Drive Train	13 months	13 months
Electronics <sup>2</sup>	13 months	13 months
Batteries	12 months	12 months

1. Battery range at maximum weight capacity and will vary due to rider weight, drive surface, terrain, and battery type.
2. Electronics warranty excludes batteries.

Literature is current at the time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

# ASSEMBLY

Your *Compass HD™* is shipped partially disassembled in order to maximize the protection of all its parts during the shipping process. Please follow the instruction below to quickly and easily assemble the power chair for your use.

## **NOTICE**

You will need only basic tools. If you do not have the required tools, or if you do not feel capable of safely assembling your power chair, please contact your local Golden Technologies, Inc. representative.

### **Main Components**

1. Frame with power module, controller cable & connectors for battery cables
2. Battery box cover
3. Seat assembly
4. Joystick
5. Batteries (Two – NF22)



**Joystick**



**Batteries**


# ASSEMBLY

## JOYSTICK INSTALLATION

1. Position the joystick on the joystick bracket. Match the holes in the joystick with the holes in the bracket. See figure 1.
2. Insert the 4 Allen head screws with lock washers through the bracket and into the holes in the joystick.
3. Use an Allen wrench to tighten the screws.
4. Connect the joystick cable to the controller cable. See figure 2.
5. Use the provided Velcro tie-wraps and under arm clip to secure the controller cable to the armrest bracket.



**Figure 1**



**Figure 2**

# ASSEMBLY

## SEAT REMOVAL/INSTALLATION



WARNING

The standard (20 x 20) seat assembly weighs about 51 pounds. Please ask for help if you do not feel capable of safely lifting that much weight. Removal of arms will significantly reduce the weight of seat.



WARNING

Be sure the seat is correctly installed and locked before operating your power chair.



WARNING

Never pull on the arms or use them as a handle to move the power chair as they could inadvertently pull out of the receiver bar causing injury.

## REMOVAL

1. Disconnect the controller cable from the joystick. “Unclip” from arm and remove Velcro tie-wraps. See figure 2 on page 13.
2. Lift the front and rear seat lock handles and slide the seat until all (4) of the seat guides are out of the receivers.

## INSTALLATION

The seat can be installed from either side of the power chair.

1. To install the seat, align the seat guides (front and rear) with the seat receivers (front and rear).
2. Slide the seat towards the center of the power chair until the seat locks securely in place. The front and rear locks must engage.
3. Reconnect the controller cable to the joystick cable.



Figure 3



Figure 4



Figure 5

Rear Seat Lock Handle  
Front Seat Lock Handle

Seat Guides  
Seat Height & Angle Adjustment

Front Seat Lock Handle

Seat Locking Handle  
Seat Receiver



Figure 6

# ASSEMBLY

## Accessing the batteries



After removing the seat, locate the (2) "Battery Box Cover Locks".



Depress both locks and lift the cover.



**Note:** On the inside of the battery box cover there is a battery connection diagram.

# ASSEMBLY

## BATTERY INSTALLATION

**NOTICE** The batteries weigh 37 pounds each. Please ask for help if you do not feel capable of safely lifting that much weight.




Figure 7




Figure 8




Figure 9



### WARNING

**Make certain that the controller power is turned off and that the chair is NOT in freewheel mode before attempting to perform assembly or disassembly.**

1. Place the batteries in the battery box. Be certain to position the terminals on each battery toward the front of the unit as shown in Figure 7. Secure the batteries in place using the hook and loop straps.
2. Use the screws, washers, and nuts provided to connect the positive (red booted) connector of each cable to the positive (+) terminal of each battery. See figure 8.
3. Use the screws, washers, and nuts provided to connect the negative (black booted) connector of each cable to the negative (-) terminal of each battery.
4. Connect the connectors of each battery cable to the corresponding connector on the main wire harness. See figure 9.
5. Cover the battery terminals with the protective caps provided. See figure 9.

# COMFORT SETTINGS

You may be spending a great deal of time on your *Compass HD™*. To provide you with the maximum seating comfort, Golden Technologies, Inc. has designed this power chair to incorporate the following adjustments for operator comfort.

## COMFORT SETTINGS

- |                    |                           |                            |
|--------------------|---------------------------|----------------------------|
| 1. Seat height     | 5. Footrest height        | 9. Joystick bracket length |
| 2. Headrest height | 6. Extending the footrest | 10. Armrest width          |
| 3. Backrest angle  | 7. Armrest angle          |                            |
| 4. Footrest angle  | 8. Armrest height         |                            |



### WARNING

**Make certain that the power to your *Compass HD™* is turned off before making any adjustments, to eliminate the risk of the joystick being accidentally bumped and activating the power chair.**

### NOTICE

If you do not feel capable of safely making these adjustments, please contact your local Golden representative.

### Seat height adjustment


The seat frame is a 4-post design that is adjustable in height and seating angle. Adjustment is made by removing ball detent pins, raising or lowering seat posts to desired height/angle and re-installing the ball detent pins.



**\*Front/Rear height difference can be 3 inches.**

# COMFORT SETTINGS

## Adjusting the seat with a front to rear angle



Loosen seat guide adjusting bolts using a 6mm Allen wrench provided in the charger bag and a pair of pliers or a 13mm open end wrench. Place the seat receivers in the lowest position and insert the ball detent pins through the towers and seat receiver legs. Leave the seat guide adjusting bolts loose. Slide the seat into the receivers, adjust seat to the desired height and angle using the (4) detent pins. Leave the seat in the seat receivers and tighten the seat guide adjusting bolts. This process is more easily accomplished with two people.




# COMFORT SETTINGS

## HEADREST HEIGHT ADJUSTMENT

1. Push and then release the clamp on the left post of the headrest while you are pulling up or pushing down on the head rest. See figure 10.
2. The headrest will “click” to a stop in one of the preset positions.
3. Repeat step 1 until the headrest is at the desired height.



**Figure 10**



**Figure 11**

## BACKREST ANGLE ADJUSTMENT


**NOTE:** The backrest can be folded down to minimize chair height during transport.

1. Pull up on the adjustment lever (see figure 11) and push against the backrest until the backrest is in the desired position.
2. Release the lever.
3. To return the backrest to the upright position, pull up on the lever and move the backrest to the desired position.

# COMFORT SETTINGS

## FOOTREST ANGLE ADJUSTMENT

1. Fold the footrest upward for easy access to the angle adjustment bolt. See figure 12.
2. Use a wrench to turn the adjustment bolt counter-clockwise to increase the footrest angle.
3. Use a wrench to turn the adjustment bolt clockwise to decrease the footrest angle.



**Figure 12**

## FOOTREST HEIGHT ADJUSTMENT

1. Remove front cover.
2. Use an Allen wrench to remove the bolts indicated in figure 13.
3. Slide the footrest to the desired height.
4. Align the bolt holes in the footrest and the footrest bracket.
5. Install and tighten the bolts.



**Figure 13**


# COMFORT SETTINGS

## EXTENDING THE FOOTREST

Refer to the pictures below.

1. Insert a 5mm Allen wrench from the top side into the (2) adjustment screws.
2. Place a 10mm wrench on the backside onto the (2) nuts and loosen the (2) footrest adjustment screws.
3. Slide the footrest plate to the desired position and tighten the adjustment screws.


**NOTE:** The maximum extended length is not at the end of the slots.



## ARMREST ANGLE ADJUSTMENT

Refer to the picture below.


1. Lift armrest up.
2. Turn jam nut counter-clockwise to loosen.
3. Turn the adjustment bolt using a 6mm Allen wrench to raise or lower the armrest-angle. Turning the bolt counter-clockwise will raise the angle. Turning the bolt clockwise will lower the angle. Turn until desired angle is reached.
4. Re-tighten jam nut by turning clockwise.
5. Repeat steps 1-4 on the other armrest.



# COMFORT SETTINGS



**Figure 14**



**Figure 15**



**Figure 16**

## **ARMREST HEIGHT ADJUSTMENT**

1. Loosen “Armrest Height Adjustment Knob” (Figure 14) until you can pull it outward.  
Note: The knob is retained in its housing with a spring.
2. While keeping the knob pulled, raise or lower the armrest to the desired position then release the knob. After releasing the knob, you may have to raise or lower the armrest until it aligns with the nearest locking position.
3. Tighten the height adjustment knob.

## **JOYSTICK BRACKET LENGTH ADJUSTMENT**

1. Loosen the setscrew that adjusts the joystick position. See Figure 15.
2. Slide the joystick bracket in or out of the joystick holder to desired position.
3. Tighten the setscrew.

## **ARMREST WIDTH ADJUSTMENT**

1. Loosen the adjustment knobs at the rear of the seat. See Figure 16.
2. Slide the arms in or out.
3. Tighten the adjustment knobs. See Figure 16.

# OPERATION



**MANDATORY!** Read and follow the owner's manual.

Your *Compass HD™* is simple to operate. However, for your safety and the safety of others, Golden Technologies, Inc. recommends that you carefully read and understand the following operating instructions. We also recommend that you practice operating your *Compass HD™* in an area free of any obstacles. Once you have gained confidence in your ability to control your power chair, you will more easily be able to operate it in normal daily conditions.

## CAUTION

Before turning on the power of your *Compass HD™*, take note of your environment and set your speed control accordingly. (See “The Speed Control Buttons” Figure 17) For indoor driving, we recommend that you select the lowest speed setting. For outdoor driving, we recommend that you select a speed setting at which you feel comfortable, safe, and in control of your power chair. Familiarize yourself with the features of your *Compass HD™* described below and follow the in-structions to safely operate your power chair.



**Figure 17**



**WARNING** The surface of the joystick can potentially get hot when it is exposed to strong sunlight for long periods.

# OPERATION

## DRIVING: Joystick Power On/Off Button



Power Off

To turn on the power to your *Compass HD™*, press the power button. If the system is healthy, the Status indicator (through the Power Button) will light up green, and the Battery Gauge will display the current battery status.

If there is a fault with the system when powering up, the status indicator will indicate a fault with a series of red flashes (see “Diagnostics”) on page 26. If the fault is one that prevents the power chair from driving, then the battery gauge will flash continuously.



Power On

To turn off the power to your *Compass HD™*, press the power button; the system will power down and the Status indicator will switch off.

## NOTICE


In the unlikely event that the power chair is in a runaway situation, press the power button to perform an EMERGENCY STOP.

## Drive Inhibit Indication

Drive inhibit mode is indicated by the battery gauge with a right-to-left light sequence.

The light sequence starts with the green LED on the right-hand side, and one-by-one, each LED will switch on then off. When the sequence completes at the left-most red LED, it begins again at the right-hand side.

The light sequence continues until the error condition has been cleared.



Drive Inhibit Light Sequence

## OONAPU

OONAPU (“**Out Of Neutral At Power Up**”) is a safety feature that prevents accidental movement of the wheelchair, either when powering up, or when the wheelchair comes out of an inhibit state.

If the LINX LE System is turned on (or comes out of an inhibit state) while the joystick is not in the center position, an OONAPU **warning** is displayed. During an OONAPU warning, the battery gauge LEDs will flash continually to alert the user, and the chair will not drive. If the joystick is returned to the center position within five seconds, the warning will clear and the wheelchair will drive normally.

However, if the joystick remains out of neutral for longer than five seconds, an OONAPU **error** will occur; the error is displayed by the Status indicator flashing red, and the chair will not drive. To clear the error, return the joystick to the neutral position and power the unit off and then on again.

# OPERATION

## The Joystick

The joystick (see figure 14 on page 21) controls the speed (up to the maximum limit set by the speed dial) and direction of your power chair.

Pushing the joystick away from the neutral (center) position will move your *Compass HD™* in the direction that the joystick is pushed.

- ◆ The farther forward or backward you push the joystick, the faster your power chair will go.

To operate your power chair, gently push the joystick in the direction in which you want to travel.

- ◆ Gentle operation of the joystick will provide you with smoother changes in speed and direction.
- ◆ Sharp or jerky operation of the joystick will result in quick and drastic changes in direction and speed.

### NOTICE

For your safety, the controller automatically sets the reverse speed, acceleration, and deceleration in proportion to the speed control setting.

### NOTICE

When you are not pushing on the joystick, or when you release the joystick, the joystick will automatically return to the neutral position, the chair will decelerate, as the electromagnetic brakes are applied, and come to a smooth stop.

## Controlling Maximum Speed



**Speed Dial**

The speed dial allows you to limit the maximum speed of the power chair (that is, the speed when the joystick is fully deflected) to suit your preference and environment.

The dial offers 10 discrete steps between the lowest speed (dial set to the left) and the highest speed (dial set to the right).



**Speed Symbol**

As a visual reminder, a speed symbol (shown left) is positioned just above the speed dial to indicate the low and high positions of the speed dial.



**Horn**

The horn button is located below the power button. Press the horn button to sound the horn. The horn will sound for as long as the horn button is pressed.

# OPERATION

## The Battery Gauge

The battery gauge comprises five different LEDs (1 x RED, 2 x AMBER, 2 x GREEN), located above the horn button. The number of LEDs lit depends on the status of the battery, as shown.

Battery Gauge



Battery Gauge

The battery gauge LEDs are also used to display charging information. Refer to “Battery Charging” on pages 30-31.

## Normal Operation

### Battery Gauge

### Battery Level



Fully charged




Consider charging the battery



Battery needs charging

## High voltage warning



A high voltage warning is indicated by all LEDs on, and the green LEDs flashing. This occurs when the battery voltage level has risen above the high voltage warning set-point.

## Low voltage warning



A low voltage warning is indicated with the left-most LED flashing. This occurs when the battery voltage level has dropped below its low voltage warning set-point. **Charge the batteries immediately.**

## Cut-off voltage



When the battery voltage drops below the battery cut-off voltage: The first (red) LED will flash on the battery gauge. The status indicator (under the power button) will display a flash code 2. The horn will sound once every 10 seconds.

# OPERATION / DIAGNOSTICS

*For Compass HD Units with the LiNX Controller*



If, when powered up, there is an error with the system, then the status indicator will flash red. The number of flashes will indicate the type of error. These are described in the table below.


## Status Indicator

Flash Code	Description	
1	Remote / joystick error	Check LiNX connections and wiring. Replace the LiNX Control Unit.
2	Network or configuration error	Check LiNX connections and wiring. Replace the LiNX Control Unit.
3	Left motor error	Check the left motor, connections and cabling.
4	Right motor error	Check the right motor, connections and cabling.
5	Left park brake error	Check the left park brake, connections and cabling.
6	Right park brake error	Check the right park brake, connections and cabling.
7	Module error (other than Remote)	Check LiNX connections and wiring. Replace the Power Module.


## NOTICE

The LiNX modules are not user serviceable. Do not attempt to repair LiNX modules.

# OPERATION



**Figure 18**



**Figure 19**

## **FREEWHEEL MODE**

- ◆ To disengage the gears and put your power chair in freewheel mode, turn the freewheel levers as shown. See figure 18.
- ◆ To reengage the gears and take your chair out of freewheel mode, turn the levers as shown. See figure 19.



**WARNING**

**Never put your power chair in freewheel mode when it is on a slope or incline of any type.**



**WARNING**

**Never put your power chair in freewheel mode while you are operating your power chair.**


## **THERMAL ROLLBACK**

Your *Compass HD™* is equipped with a safety system. A microprocessor monitors the operating temperatures of the controller. In the event of excessive heat occurring in the controller, the controller will decrease the speed of your chair. This is done to reduce the load on the electrical system and allow the components to dissipate heat. The controller will automatically set the chair's speed back to full normal when the operating temperature returns to normal levels.

## **THE MAIN CIRCUIT BREAKER**

The main circuit breaker is another safety feature incorporated into your *Compass HD™*. This device monitors the amount of current being drawn from the batteries when the chair is in use. When the motors are heavily strained and too high a current draw is being placed on the batteries, the main circuit breaker will trip and bring your chair to a stop.

# OPERATION



**Figure 20**

## **Reset Button**

If the main circuit breaker trips, please wait for approximately 1 minute and then push the reset button to reset the main circuit breaker. See figure 20.

**Note:** Usually, the thermal rollback feature is more sensitive than the main circuit breaker. We recommend that you turn off the power to the chair and wait for five minutes when the chair suddenly loses speed or power. Doing so will allow the overheated electrical components to cool to their normal operating temperatures.

## **SWING-AWAY DESK ARM**

Your *Compass HD™* is equipped with a swing-away desk arm that provides you with the convenience of positioning your chair close to a work surface such as a table or desk.

To swing the desk arm simply push on the side of the joystick until the joystick swings to the retracted position. See figures 21, 22, and 23.



**Figure 21**



**Figure 22**



**Figure 23**

# BATTERY CHARGING



## WARNING

**You must use only the charger that is supplied with your *Compass HD™*. The use of any other charger on this power chair will void the warranty. Using unauthorized chargers may also result in severe damage to the batteries and/ or damage to the chair. Using the wrong charger may also be a hazard.**

## Operational Range

Depending on the use, the terrain, and the driving conditions, the batteries of your *Compass HD™* will provide approximately twelve miles of operation.



**WARNING!** Corrosive chemicals contained inside battery.  
Explosive battery conditions exist.

**WARNING!** Charger surface can become HOT.  
Electrical hazard!

## BATTERIES AND CHARGING

Battery maintenance is the most important part of maintaining your power chair. Keeping the batteries fully charged helps to extend battery life. Use the following guidelines to help keep your batteries in optimum condition.

**IMPORTANT!** New batteries **MUST** be fully charged prior to the initial use of your power chair.

**Charge new batteries for 12 hours even if the battery gauge already shows a full charge. This is to condition your new batteries for maximum performance.**

- For daily use, keep batteries fully charged. We recommend that you plug in the off-board charger after each use and charge 6-8 hours.

**“For best results please charge your batteries overnight after every use.”**

- If you are not going to use the power chair for more than a week, **fully charge** the batteries and then disconnect them from the power chair. Refer to page 16.

## NOTICE

### **Charging Guidelines Checklist to Maximize Battery Life**

- ✓ Use only the automatic off-board charger supplied for all routine charging.
- ✓ Never use an automotive or wet type charger.
- ✓ Avoid deep discharges and never drain the batteries completely.
- ✓ Do not leave batteries in a low state of charge for extended periods. Charge a discharged battery as soon as possible.
- ✓ Fully recharge batteries regularly.
- ✓ Always store batteries fully charged.
- ✓ Check stored batteries once a month and recharge as necessary.

# BATTERY CHARGING

## CHARGING THE BATTERIES

Your *Compass HD*'s charging system is designed for your safety and for your convenience. Follow the steps below to recharge the batteries.

1. Position your *Compass HD*<sup>TM</sup> close to a standard wall electrical outlet.
2. Turn off the power on the joystick.




**WARNING** The charger should not be plugged into the electrical outlet or powered up at this point in the process. Please check before proceeding.

4. Insert the matching charger plug into the joystick's charging socket that is located at the front bottom edge of the joystick. See figure 24.

**IMPORTANT: Failure to align pins before inserting the plug in the correct orientation may damage the joystick and render your wheelchair non-operational.**

5. Insert the plug at the other end of the charger power cord into a standard electrical wall outlet.
6. Disconnect the charger power cord from the wall outlet and from the joystick charging socket or the battery pack for (alternative charging method) when the batteries are fully charged.

**NOTICE** The batteries should be fully charged in approximately eight hours. A full charge is indicated when the battery charger LED is green and the 5 battery gauge LEDs (1 x RED, 2 x AMBER, 2 x GREEN) on the joystick are illuminated.




**Figure 24**

## Battery Charging

The battery gauge will indicate the system is being charged by cycling between a left-to-right light sequence, and then display the approximate battery charge state at the end of the sequence.

Driving is prevented (inhibited) while the system is being charged.

The system does not have to be powered up when charging the battery, however, if it is not powered up, then the battery gauge will not display the charging state/light sequence.



# CARE AND MAINTENANCE

## ROUTINE MAINTENANCE

Your power chair requires a minimal amount of care and maintenance. If you do not feel confident in your ability to perform the maintenance listed below, you can schedule an inspection and maintenance with your provider.

- ✓ Make sure your batteries are fully charged daily. See “Battery Charging”.
- ✓ While turned off, check the joystick control daily for damage, and make sure it returns to the neutral position. If any defects are observed, contact your provider. **Do not attempt to repair it yourself!**
- ✓ Inspect the seating system, armrests, footplate, and front riggings daily for loose hardware, excessive wear, or damage. If any defects are observed, contact your provider.
- ✓ Regularly inspect your power chair’s tires for signs of wear.
- ✓ Regularly inspect all exposed harnessing and connections for wear, damage, and corrosion. Have your provider repair or replace any damaged harness or connector.
- ✓ Keep all the electrical components (charger, joystick control, and electronics) free from moisture. If any of these items do become exposed to moisture, let them dry thoroughly before attempting to operate your power chair again.
- ✓ Periodically check all fasteners for tightness, wear, and corrosion. If any of these are present, turn off the power chair and contact your provider.



WARNING

**Be careful of pinch points and sharp metal parts (wear gloves).**

- ✓ All wheel bearings are permanently lubricated and sealed. No additional lubrication is required.

We recommend that you check the following:

**Tire Tread:** Regularly visually inspect the tire tread. If the remaining tread is less than 1/32 of an inch, have your local Golden representative replace the worn tires.

**Joystick and Controller:** Protect the joystick and controller from adverse weather conditions. Moisture will damage the controller and void your power chair’s warranty.

**Cleaning:** Use only a water-dampened soft cloth to clean the ABS shroud. The shroud has a protective coat that is very easy to maintain.



WARNING

**Do not use free-flowing water to clean your Compass HD™. Water and extreme temperatures are the main elements that can adversely affect your power chair and its performance.**

### **Water, Rain, Sleet, and Snow**

Water, in any form, will cause electronic malfunction or corrosion of the electrical components, connections, and the chair frame.

### **Temperature**

- ◆ At extremely low temperatures, the batteries of your power chair may freeze, preventing your *Compass HD™* from operating.
- ◆ At extremely high temperatures, your power chair may operate at slower speeds due to the controller’s thermal rollback feature that is designed to prevent damage to the motors and to the other electrical components of the chair.

# DISASSEMBLY



## WARNING

**Make certain that the controller power is turned off and that the chair is NOT in freewheel mode before attempting to perform disassembly.**

### SEAT REMOVAL

1. Disconnect the joystick cable from the joystick. See figure 2 on page 14.
2. See page 15 to complete the removal.



## WARNING

**The standard (20 x 20) seat assembly weighs about 51 pounds. Please ask for help if you do not feel capable of safely lifting that much weight. Removal of arms will significantly reduce the weight of seat.**

### BATTERY REMOVAL for Transportation

After removing the seat and battery box cover:

1. Disconnect the battery connectors and unfasten battery straps. See figures 7 and 9 on page 16.
2. Use the battery handles to lift the batteries out, one at a time.



## WARNING

**The batteries weigh 37 pounds each. Please get help if you do not feel capable of safely lifting that much weight.**

### For Replacement

After removing the seat and battery box cover:

1. See pages 10 and 11.
2. Disconnect the battery connectors and unfasten battery straps.
3. Use the battery handles to lift the batteries out, one at a time.
4. Remove the protective caps from the battery terminals.
5. Disconnect the positive (+) battery cables from the battery terminals, one at a time.
6. Disconnect the negative (-) battery cables, one at a time.
7. Install new batteries. See BATTERY INSTALLATION on page 16 or see the instructions on the inside of battery box lid and follow the steps.

# WARRANTY



## GOLDEN

### Consumer Limited Warranty *Power Wheelchairs*

Golden Technologies has provided information throughout the owner's manual in regards to proper usage and care of your power wheelchair. Failure to adhere to this information will void the warranty of the power wheelchair, its parts and the batteries supplied by Golden Technologies

The following warranty is granted only to the INITIAL Consumer who has purchased our product and commences on the Date of Purchase by the Consumer from an Authorized Golden Technologies' Provider. The Consumer must have a warranty registration card filed with Golden Technologies to receive service.

The warranty covers the following power wheelchair models unless otherwise noted:

- ◆ LiteRider Envy GP162
- ◆ Compass Sport GP605
- ◆ Alante Sport' GP208
- ◆ Compass HD GP620

#### *Five-Year Limited Warranty:*

For the period of five years from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will repair or replace at our option any of the following **structural frame components** found to be defective by an authorized Golden Technologies representative:

- ◆ Steel Frame\*
- ◆ Frame Welds\*
- ◆ Seat Post\*

*\*The GP605 is warranted for LIMITED LIFETIME*

#### *Thirteen Month Limited Warranty*

For the period of thirteen months from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will repair or replace at our option any of the following **structural frame components** found to be defective by an authorized Golden Technologies representative:

- ◆ Caster Beam
- ◆ Caster Forks
- ◆ Footrest
- ◆ Anti-Tip Forks
- ◆ Armrests

For the period of thirteen months from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will repair or replace at our option any of the following **drive train components** found to be defective by an authorized Golden Technologies representative:

- ◆ Motor/Gearbox Assembly
- ◆ Brake Assembly (electrical function only)

**Note: An increase in operating noise of the motor/gearbox DOES NOT constitute a defect. With normal wear and tear, operating noise is expected to increase.**

For the period of thirteen months from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will repair or replace at our option any of the following **electronic assemblies** found to be defective by an authorized Golden Technologies representative:

- ◆ Electronic Controller
- ◆ Battery Charger
- ◆ Wire Harness(es) (excludes boot covers)

**Note: While charging batteries, it is normal for the battery charger to heat up. Heat coming from the battery charger DOES NOT constitute a defect.**

For the period of thirteen months from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will repair or replace at our option any of the following **parts** found to be defective by an authorized Golden Technologies representative:

- ◆ Bearings
- ◆ Rubber Components
- ◆ Plastic Components
- ◆ Bushings
- (excludes tires and all wear parts)
- (excludes body & all wear parts)
- ◆ Seat swivel mechanism

#### *Twelve Month Limited Warranty*

For the period of twelve months from the date of purchase, in the event of defective materials or workmanship, Golden Technologies will replace the **batteries supplied by Golden Technologies** if found to be defective by an authorized Golden Technologies representative.

# WARRANTY



## Consumer Limited Warranty *Power Wheelchair (continued)*

Items not mentioned or listed are covered at the discretion of Golden Technologies. Any accessories, standard or optional, supplied by Golden Technologies, are covered for a period of one year from the date of purchase with the particular unit.

### ***Warranty Exclusions:***

Golden Technologies does not provide warranty on any of the following items which may require replacement due to the normal wear and tear of day to day usage:

- ◆ Tires and Tubes
- ◆ ABS Plastic Shrouds
- ◆ Armrest Pads
- ◆ Motor Brushes
- ◆ Bulbs/Fuses
- ◆ Battery Cases
- ◆ Brake Pads
- ◆ Upholstery

This warranty also excludes the following:

- ◆ Loss or theft of components
- ◆ Damage caused by:
  - Battery fluid spillage or leakage.
  - Abuse, misuse, accident or negligence.
  - Improper operation, maintenance or storage.
  - Commercial use or use other than normal.
  - Repairs and/or modifications made to any part without the specific consent of Golden Technologies.
  - Exceeding the specified weight capacity of the unit.
  - Accessories other than those supplied or approved by Golden Technologies.
  - Failure to adhere to the product instructions.
  - Acts of Nature, such as lightning strikes, earthquakes, etc....
  - Circumstances beyond the control of Golden Technologies.
- ◆ **ANY PARTS ALTERED OR REPAIRED BY UNAUTHORIZED PERSON.**

### ***Warranty Service:***

Warranty Service must be performed by an authorized Golden Technologies representative. Golden Technologies reserves the right to replace warranted part(s) with refurbished or new part(s) at our discretion. All labor charges, service calls, transportation costs or any other charge(s) associated with the installation of any warranted part(s) are the responsibility of the consumer. Warranty is for the replacement of the parts only and does not include freight for the replacement parts. Consumers are not to return any item(s) to Golden Technologies without prior written authorization. Any damages incurred while warranted part(s) are in transport are the sole responsibility of the consumer.


### **There is no other express warranty.**

Implied warranties, including those of merchantability and fitness for a particular purpose are excluded. Liabilities for consequential damages are excluded.

This warranty gives you specific rights and you may also have other rights which may vary from state to state.



Fold Here First



Place  
Stamp  
Here

Golden Technologies  
401 Bridge St.  
Old Forge, PA 18518



Fold Here Second



401 Bridge Street, Old Forge PA 18518  
Tel: (800) 624-6374 Fax: (800) 628-5165  
[www.goldentech.com](http://www.goldentech.com)