

AL435 3-AXIS

THE ULTIMATE HOIST LIFT SOLUTION



PREMIUM 3-AXIS INSIDE LIFT

The AL435 powered-boom maneuvers even large power wheelchairs or scooters into small vehicle openings. More versatile and flexible, for use inside a van, SUVs and pick-up trucks. The AL435 provides fully powered lifting, boom extension, and rotation for the easiest and most effortless loading. The power extend function allows for precise positioning while rotating into a tight cargo area.

- Power lifting, extension, and rotation
- Simplified base fits in a smaller footprint
- Passenger or driver side installation
- Weather-resistant hand control
- Universal base adapter included



PENDANT CONTROLLER

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics



We wanted a lift that would fit in our minivan and be easy for us to operate when loading and unloading my husband's scooter. The AL435 was our first choice, we are able to retain use of the second row of seats, which is nice when we have our grandkids.

AL435 3-AXIS PRODUCT SPECIFICATIONS

LIFTING CAPACITY 400 lb

INSTALLED WEIGHT 76 lb

HORIZONTAL REACH 21" - 32"

ADJUSTABLE HEIGHT 32" - 43"

WARRANTY 3-year

SECUREMENT Docking device *(not included)*

INCLUDED

- Folding arm
- Universal base adapter

ACCOMMODATES

- Manual wheelchairs
- Power wheelchairs
- Scooters

POWERED

- Power lifting
- Power rotation
- Power extension

AVAILABLE OPTION(S)

- Battery pack
- Custom adapters
- Custom docking devices
- Wireless remote