

# AL825

## PICKUP 225 SPEED LIFT



### PICKUP 225 SPEED LIFT

Fits regular or extended cab trucks for manual wheelchairs, small power wheelchairs, and scooters. Speed model is two times faster than heavy-duty models and ideal for manual wheelchairs.



#### PENDANT CONTROLLER

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics

#### FEATURES:

- Delivers mobility device close to driver or passenger side door
- Powered lifting and rotation
- Offset base mounting fits tightly to sidewall and fender walls
- Driver or passenger-side mounting bases available
- Dashboard-mounted ON/OFF switch
- Dual-strapped hook prevents mobility device from spinning during loading/unloading
- Pendant controller included (*wireless option available*)

### AL825 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	210 lb	<b>INSTALLED WEIGHT</b>	116 lb
<b>ADJUSTABLE LENGTH</b>	28.75" - 46.5"	<b>ADJUSTABLE HEIGHT</b>	34" - 41"
<b>SECUREMENT</b>	Docking device ( <i>not included</i> )	<b>INCLUDED ACCESSORY</b>	Universal base adapter
<b>APPLICATION</b>	<ul style="list-style-type: none"><li>- Manual wheelchairs</li><li>- Small power wheelchairs</li><li>- Scooters</li></ul>	<b>POWERED</b>	<ul style="list-style-type: none"><li>- Power lifting</li><li>- Power rotation</li></ul>
<b>AVAILABLE OPTIONS</b>	<ul style="list-style-type: none"><li>- Battery pack</li><li>- Custom adapters</li><li>- Wireless remote</li></ul>	<b>WARRANTY</b>	<ul style="list-style-type: none"><li>- 3-year</li><li>- Transferable</li></ul>