



FEATURES & BENEFITS

- Folds for easy storage
- **10" to 12" occupied interior space**
- 1,000 lb load capacity
- All-Weather Traction
- Mounts to floor of many standard vans (Check for compatibility)



Manual Folding Van Ramps

Installs into many standard factory vans. A more cost-effective alternative to a fully converted accessible van. Operated by raising and lowering the ramp from a folded position to the ground.

- Ramp is raised to an upright position inside the van door.
- Secured with included security strap.
- Ships as a complete unit with hardware.
- Occupies approximately 10" of depth.

Power Folding Van Ramps

An alternative to a lowered floor conversion. Operated with the push of a button. Ramp deploys from a folded position inside the van down to the ground.

- Built-in safety feature: Ramp only moves when button is being pressed.
- Reliable direct-drive actuator motor.
- Ships as a complete unit with hardware.
- Pendant cord or wireless remote control options.
- Occupies approximately 12" depth in folded position.

Available ramp widths:
26"
30"
36"

Also, See The Flyer
For Our
Assisted Van Ramps!

Available ramp widths:
26"
30"

REAR INSTALLATION



SIDE INSTALLATION






CHECK COMPATIBILITY WITH VAN PRIOR TO ORDERING

1) Measure Door Width and Height 2) Measure Wheelbase of MoDevice 3) Refer to Chart Below

MINI VAN MANUAL RAMPS			
		DEALER	RETAIL
MF-26	6' x 26" Manual Folding Ramp 40" door height		
MF-26.7	7' x 26" Manual Folding Ramp 46" door height		
Requires 28" minimum side door opening, maximum chair wheelbase is 25"			
MF-30.6	6' x 30" Manual Folding Ramp 40" door height		
MF-30.7	7' x 30" Manual Folding Ramp 46" door height		
Requires 32" minimum side door opening, maximum chair wheelbase is 29"			
FULL SIZE VAN MANUAL RAMPS			
MF-30.7	7' x 30" Manual Folding Ramp 46" door height		
MF-30	8' x 30" Manual Folding Ramp 52" door height		
MF-30.9	9' x 30" Manual Folding Ramp 58" door height		
Requires 32" minimum side door opening			
MF-36.7	7' x 36" Manual Folding Ramp 46" door height		
MF-36	8' x 36" Manual Folding Ramp 52" door height		
MF-36.9	9' x 36" Manual Folding Ramp 58" door height		
Requires 38" minimum side door opening, maximum chair wheelbase is 35"			
POWERED MINI VAN RAMPS			
		DEALER	RETAIL
MF1	6' x 26" Power Mini Fold Pendant Cord		
MF1-R	6' x 26" Power Mini Fold Wireless Remote		
Door opening: minimum 29" x 42", maximum chair wheelbase is 25"			
POWERED FULL SIZE VAN RAMPS			
AF1	7' x 30" Power Auto Fold Pendant Cord		
AF1-R	7' x 30" Power Auto Fold Wireless Remote		
Door opening: minimum 33" x 48", maximum chair wheelbase is 29"			



**VAN RAMPS ARE SHIPPED FREIGHT:
QUOTE REQUIRED**

Roll-A-Ramp® Return Policy

*20% Restocking Fee

*Return within 90 days from purchase

*Outbound shipping is not refunded

*Customer responsible for ramp return

Refund issued upon return of items in acceptable condition.

Need more weight capacity? Roll-A-Ramp HD can hold up to 1,500 pounds!