



Redefining standards in care



A new adjustable electric bed

Safe, Versatile, Efficient



Yes, you can.®

Promoting comfort and rest

Introducing the new generation of CS Series beds from Invacare.

Our new CS600™ adjustable electric bed promotes comfort and rest for patients at home or in long-term care environments.

With a focus on safety and on delivering our highest quality standards to ensure all of our beds support patients as they rest.



Redefining standards in care

Introducing the new CS600™ adjustable electric bed, the next generation of our popular CS Series of beds.

Built from a design foundation of exceptional comfort, safety and versatility, the CS600™ design is focused on improving standards in patient care.

Safety first. No compromise.

When it comes to patient and caregiver safety there is no compromise. The CS600™ bed meets the highest international safety and performance standard.*

*Medical Bed Equipment Standard for safety, quality and performance - IEC 60601-2-52.



1



Safe

Reassurance in quality is at the core of what we do. The CS600™ is tested to the highest safety and performance standard* for peace of mind.

2



Versatile


The versatile width design of the CS600™ and up to 600lb safe working load meets a varied profile of patients and care environments.

3




Efficient

Tool-free assembly reduces setup time and and complex installation. Slat deck design and wipe down electronics make cleaning a breeze.



Height adjustable
The CS600™ lowers to the floor. A lower bed height may help reduce the possibility of patient falls. For caregiver convenience, the bed also reaches a height of 30".




Repositioning the bed to a cardiac upright position helps to improve circulation and general healing.



Easy to clean slat deck design helps to promote improvements in infection control.



Choice of bed ends
Available in 3 finishes



Safe working load 600 lbs



Backlit hand control for improved visibility




One touch floor lock
Ease and efficiency for staff

Accessories available:



Length extension between 80" to 84"



Width extension set to 36", 39" or 42"



Trapeze secures in 3 positions to accommodate patient limitations. 200 lb safe working load



Trolley for easy transportation



Battery backup for operation in case of power outage



Underbed light for improved floor visibility

Technical data ▼

					
Bed Deck Length ▼	Bed Deck Width ▼	Overall Length (incl. head/foot boards) ▼	Maximum Safe Working Load ▼	Minimum Bed Deck Height ▼	Maximum Bed Deck Height ▼
80"	CS600-36 36" CS600-EXP 36", 39", 42"	88"	600 lbs	7"	30"
					
Maximum Back Angle ▼	Maximum Knee Angle ▼	Maximum Reverse Trendelenburg ▼	Length Extension ▼	Width Extension ▼	ACP (Attendant Control Panel) ▼
60°± 5°	30°± 5°	17°± 5°	80" - 84"	CS600-36 w/accessory 42"	Optional

Bed models

Standard Width: 36" W by 70" L
Item #: CS600-36

Bed end colors



Solar Oak



African Walnut



Biltmore Cherry

Frame colors



Light Gray

Please note: colours may vary slightly from those displayed above.

Warranties

- ▶ Actuator - 2 years
- ▶ Control box - 1 year
- ▶ Hand control - 1 year
- ▶ Frame - 2 years

Reference user manual for additional warranty details.

Tel: (800) 333-6900
www.invacare.com

Invacare Corporation
Worldwide Headquarters
1 Invacare Way
Elyria, OH 44035-4190



Yes, you can.®