

AL010/AL050

MICRO OUTSIDE LIFT



COMPACT SOLUTION FOR TRAVEL SCOOTERS AND WHEELCHAIRS

Specifically designed to carry compact three (3) and four (4) wheel travel scooters, manual wheelchairs, and small power wheelchairs.

AL010 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	135 lb
ACCOMMODATES	Scoters
INSTALLED WEIGHT	71 lb
PLATFORM SIZE	23.5" x 33" adjustable wheel cradle extends up to 42.25"
HITCH CLASS	I or greater
POWERED	- Manual fold - Power lifting
WARRANTY	3-year
INCLUDED	- AL107 hitch post - AL126 class II adapter - License plate mount
SECUREMENT	- 4 manual straps - Rear wheel cradle
AVAILABLE OPTION(S)	- Battery pack - Cover - License plate riser - Swing-away

AL050 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	135 lb
ACCOMMODATES	Manual wheelchairs Small power wheelchairs
INSTALLED WEIGHT	67 lb
PLATFORM SIZE	23.5" x 33"
HITCH CLASS	I or greater
POWERED	- Manual fold - Power lifting
WARRANTY	3-year
INCLUDED	- AL107 hitch post - AL126 class II adapter - License plate mount
SECUREMENT	4 manual straps
AVAILABLE OPTION(S)	- Battery pack - Cover - License plate riser - Swing-away