

AL030 POWER TOTE



IDEAL FOR MANUAL WHEELCHAIR TRANSPORT

Lightweight, powered lift designed for transporting manual wheelchairs. Features a 100 lb capacity with simple, one-switch operation. Adjustable hold-down arm safely secures the chair. Lift folds when not in use.

One-switch operation



Adjustable hold-down arm

Wheel cradle



The power tote has made getting around so much easier. Instead of struggling to get my wheelchair into the back seat, now it can simply roll onto the lift, power into position, and away I go. Saves room too!

AL030 PRODUCT SPECIFICATIONS

| | | | |
|----------------------------|--|----------------------|--|
| LIFTING CAPACITY | 100 lb | ACCOMMODATES | Manual folding wheelchairs |
| INSTALLED WEIGHT | 49 lb | PLATFORM SIZE | 14" x 30" |
| HITCH CLASS | I or greater | POWERED | Power lifting and manual fold |
| WARRANTY | 3-year | INCLUDED | Ball mount |
| AVAILABLE OPTION(S) | <ul style="list-style-type: none"> - Battery pack - Swing-away - Wheelchair cover | SECUREMENT | <ul style="list-style-type: none"> - Manual hold-down arm - Wheel cradle |