

Economy Inside Lift AL055



Economy Inside Lift

Super compact and convenient! The affordable AL055 trunk lift lifts and secures micro / travel power chairs and scooters into a cargo area or trunk.

Features & Benefits

- Lightweight and compact
- Heaviest component is just 23 lbs.
- Small enough to mount in trunk or cargo area
- Easy to disassemble for storage
- Power lifting with manual rotation
- No scooter or power chair modifications required
- Remote hand control included
- DD17 included; upgrades available



Specifications	AL055
Accommodates	Micro / travel power chairs & scooters
Power	Power lifting; manual rotation
Lifting Capacity	200 lbs. (91 kg)
Installed Weight	40 lbs. (18 kg)
Securement	Docking device
Accessories Included	Docking device 17
Options	Battery pack
Warranty	3-year transferable



Remote hand control included



Optional AL205 Battery Pack avoids hard-wiring into vehicle or serves as emergency backup.



800-833-0478 | www.harmar.com

