

AL160

PROFILE SCOOTER LIFT



SEE THE VIEW, NOT THE LIFT

The open platform design of the AL160 creates a simple, streamlined silhouette. Its understated presence on a vehicle is stylish and ideal for compact cars.

Harmar's AL160 Profile Scooter Lift safely and securely transports power scooters with a wheelbase between 30" and 42" – no scooter modifications needed. The automatic hold-down foot secures your scooter automatically with the push of a single button. Even larger power scooters can be transported with this profile lift – up to 350 lb.



The AL160 profile lift is ideal for my scooter. Its open platform frame allows for great visibility when driving without my scooter and has plenty of capacity for when I need to take it along.



AL160 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb	ACCOMMODATES	Scooters up to 42" wheelbase
INSTALLED WEIGHT	74 lb	PLATFORM SIZE	25" x 25" adjustable cradle extends up to 33"
HITCH CLASS	I or greater	POWERED	Fully powered
WARRANTY	3-year	INCLUDED	AL123 class III hitch adaptor, AL107 hitch post
AVAILABLE OPTION(S)	- Battery pack - Swing-away - Scooter cover	SECUREMENT	- Hold-down foot - Rear wheel cradle