AL300 / AL300HD FUSION LIFT





	* * =		
		∭ <u>*</u> * * *	* * *
* * *		* * * *	* * * *

CONVENIENCE AND VERSATILITY IN ONE

The AL300 features Harmar's Stars-N-Stripes aluminum deck design, dedicated to our veterans and commemorating Harmar's heritage as an American company.

Compatible with most scooters and power wheelchairs, it safely secures the mobility device with four, easy-to-operate, Q'Straint retractors. When empty, the platform automatically folds up and out of the way with the flip of a switch.

AL300 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb	
ACCOMMODATES	Power wheelchairsScooters	
INSTALLED WEIGHT	80 lb	
PLATFORM SIZE	26" x 39" adjustable cradle extends up to 46"	
HITCH CLASS	II or greater	
POWERED	Fully powered	
WARRANTY	3-year	
INCLUDED	 AL107 hitch post AL123 class III hitch adapter License plate mount 	
SECUREMENT	 4 self-tensioning Q'Straint retractors Rear wheel cradle 	
AVAILABLE OPTIONS	 Battery Pack Cover License plater riser Swing-away 	

AL300HD PRODUCT SPECIFICATIONS

LIFTING CAPACITY	400 lb	
ACCOMMODATES	Power wheelchairsScooters	
INSTALLED WEIGHT	82 lb	
PLATFORM SIZE	26" x 39" adjustable cradle extends up to 46"	
HITCH CLASS	III or greater	
POWERED	Fully powered	
WARRANTY	3-year	
INCLUDED	 AL107 hitch post AL123HD class III heavy duty hitch adapter License plate mount 	
SECUREMENT	 4 Self-tensioning Q'Straint retractors Rear wheel cradle 	
AVAILABLE OPTIONS	 Battery pack Cover License plate riser Swing-away 	



26SEP2020 | MKT-000029 REV A

