

AL425HD

STRONG ENOUGH TO GET THE JOB DONE



HIGH CAPACITY 450 LB HOIST LIFT

This 450 lb capacity lift solves the challenges created by heavier, more technologically advanced power chairs. Harmar's heavy duty construction and stronger components still allow for a light install weight of 70 lb. Strong enough to get the job done, while still not adding significant weight to your vehicle.

AL425HD BENEFITS:

- One person operation
- Powered lifting and rotation functions
- Built with offset post that allows for easy curbside loading
- Ideal for pick-up trucks and full size vans
- Passenger or driver-side compatibility
- Weather resistant hand controls or an optional wireless control



PENDANT CONTROLLER

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics



No need to continue to use ramps that take up space and require more set up time. The AL425HD does the heavy lifting for me.

AL425HD PRODUCT SPECIFICATIONS

LIFTING CAPACITY 450 lb

INSTALLED WEIGHT 70 lb

ARM HEIGHT 31.5" - 46.5"

ADJUSTABLE LENGTH 26.25" - 33.5"

POWERED

- Power lifting
- Power rotation

SECUREMENT Docking device *(not included)*

ACCOMMODATES

- Manual wheelchairs
- Power wheelchairs
- Scooters

AVAILABLE OPTIONS

- Custom adapters
- Custom docking devices
- Folding arm
- Wireless remote

INCLUDED ACCESSORY Universal base adapter

WARRANTY

- 3-year
- Transferable