

AL580 / AL580XL

NEXT GENERATION POWER WHEELCHAIR LIFT



LIGHTEST SOLUTION ON THE MARKET

The AL580 and AL580XL are compact, durable outside vehicle platform lifts that safely and securely transports power chairs.

The AL580 series is easy to operate and built to last for years of trouble-free use. Designed to meet the needs of consumers who need a quick, automatic securing option without making modifications to their power wheelchair, the stylish and compact detailing of this lift makes it a must have for power wheelchair owners.

AL580 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	Mid-wheel drive power wheel-chairs with 10" drive wheels
INSTALLED WEIGHT	63 lb
PLATFORM SIZE	16.5" x 27"
HITCH CLASS	II or greater
POWERED	Fully powered
WARRANTY	3-year
INCLUDED	<ul style="list-style-type: none">- AL107 hitch post- AL123 class III hitch adapter- License plate mount
SECUREMENT	<ul style="list-style-type: none">- Padded hold-down arm- Wheel cradle
AVAILABLE OPTION(S)	<ul style="list-style-type: none">- Battery pack- Cover- License plate riser- Swing-away

AL580XL PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	Mid-wheel drive power wheel-chairs with 14" drive wheels
INSTALLED WEIGHT	63 lb
PLATFORM SIZE	17.5" x 27"
HITCH CLASS	II or greater
POWERED	Fully powered
WARRANTY	3-year
INCLUDED	<ul style="list-style-type: none">- AL107 hitch post- AL123 class III hitch adapter- License plate mount
SECUREMENT	<ul style="list-style-type: none">- Padded hold-down arm- Wheel cradle
AVAILABLE OPTION(S)	<ul style="list-style-type: none">- Battery pack- Cover- License plate riser- Swing-away