

AL815CC + AL825CC

CREW CAB LIFTS



LIGHTWEIGHT AND ECONOMICAL

Harmar's Crew Cab lifts offer a quality choice in capacities to handle most mobility devices. Each offers powered lifting and rotation with a boom that adjusts in both height and length for more flexibility in installation and use.



PENDANT CONTROLLER

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics

FEATURES:

- Delivers chair close to driver or passenger side door
- Powered lifting and rotation
- Only two motors to control with removable hand control
- In-cab remote control extension included

AL815CC PRODUCT SPECIFICATIONS

LIFTING CAPACITY	135 lb
ACCOMMODATES	<ul style="list-style-type: none">- Manual wheelchairs- Small power wheelchairs- Scooters
INSTALLED WEIGHT	121 lb
ADJUSTABLE LENGTH	49.25" - 67"
ADJUSTABLE HEIGHT	34" - 41"
POWERED	<ul style="list-style-type: none">- Power lifting- Power rotation
WARRANTY	<ul style="list-style-type: none">- 3-year- Transferable
INCLUDED	Universal base adapter
SECUREMENT	Docking device <i>(not included)</i>
AVAILABLE OPTIONS	<ul style="list-style-type: none">- Battery pack- Custom adapters- Wireless remote

AL825CC PRODUCT SPECIFICATIONS

LIFTING CAPACITY	225 lb
ACCOMMODATES	<ul style="list-style-type: none">- Manual wheelchairs- Small power wheelchairs- Scooters
INSTALLED WEIGHT	121 lb
ADJUSTABLE LENGTH	49.25" - 67"
ADJUSTABLE HEIGHT	34" - 41"
POWERED	<ul style="list-style-type: none">- Power lifting- Power rotation
WARRANTY	<ul style="list-style-type: none">- 3-year- Transferable
INCLUDED	Universal base adapter
SECUREMENT	Docking device <i>(not included)</i>
AVAILABLE OPTIONS	<ul style="list-style-type: none">- Battery pack- Custom adapters- Wireless remote