WER EZ 2

ITEM# F-PEZPL2





ABOUT

This lift is known for its versatility and affordability, and is a popular choice for residential pools for all those reasons and more. Featuring battery powered lifting capabilities, this model also rotates manually requiring assistance from others for operation.



FEATURES:

- 400 lb (181 kg) weight capacity
- 360 degree manual rotation
- Comes with durable mesh sling*
- Adjustable stainless steel chains
- Stainless steel construction
- Durable powder coat finish
- Limited five (5) year warranty**
- Ships small parcel, some assembly required
- Made in USA

*Sling seat available in different sizes. Seat belt options also

** One year warranty on sling and hydraulic actuator.

PROUDLY SOLD BY:





STANDARD COLOR OPTIONS:



This lift does not include plastic options. Hard seat option comes in standard blue.





This lift does not include powder coat color options.

POWER EZ 2

APPLICATIONS:

This lift works on in ground and partially raised pools with clearances up to 15" (38 cm). This lift can also be anchored to pool decks made of concrete, pavers or wood. Choose the appropriate anchor system for your needs. Anchors sold separately. To see if the Power EZ 2 will work for you, visit our website and download the application guide, or contact our friendly sales team for assistance.

Standard Anchor *F-033EZA*



Paver Anchor *F-30220*



Surface Mount Anchor F-035EZSM



Wood Deck Anchor *F-30600*



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS & ACCESSORIES:

- Anchor options for concrete, paver, and wood deck applications
- Cover to protect your investment from the elements*
- Four-point hanger bar for added safety
- Hard seat option
- Solar charger to maintain and extend battery life
- Various sling styles and sizes to choose from

*Color and material options available

SCAN ME



Use your mobile device to scan the code and view the product webpage.













