

Model	SC134
Weight Capacity	500 lbs.
Maximum Speed <sup>2</sup>	Up to 9.3 mph
Ground Clearance <sup>4</sup>	3.5" to motor, 4.4" at mid-deck
Turning Radius <sup>4</sup>	82.28"
Overall Length <sup>4</sup>	55.9"
Overall Width <sup>4</sup>	27.05"
Seat-to-Ground Height	24.2" - 25.6"
Seat-to-Deck Height	16.5" - 18.4"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Full suspension
Range Per Charge <sup>2,3</sup>	Up to 23.4 miles
Total Weight	201.2 lbs. heaviest piece 253.8 lbs. with seat 298.8 lbs. with batteries 351.4 lbs. with seat and batteries
Heaviest Piece When Disassembled <sup>4</sup>	N/A (Unit does not disassemble)
Standard Seating <sup>3</sup>	Type: High back Weight: 52.6 lbs. Material: Black vinyl Dimensions: Width: 18", Depth: 18"
Drive System	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>5</sup>	Type: (2) 12V, sealed lead-acid, non-spillable deep-cycle Size: 75Ah Weight: 48.8 lbs each
Battery Charger	8 amp, off-board
Warranty	Frame: Lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6 months

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC





DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

1. Pride<sup>®</sup> FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of +10%, -5%.
3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4. Tested in accordance with ANSI/RESNA, WC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
5. Battery weight may vary based on battery manufacturer.

0005410-080423

0005410-082823

**PX4**  
MORE ROOM  
TO ROAM

FDA Class II Medical Device<sup>1</sup>





**9.3 MPH**  
CRUISING SPEED

**FULL**  
SUSPENSION

**500 LBS**  
WEIGHT CAPACITY

**23.4 MILE RANGE**

**LIVE YOUR BEST<sup>®</sup>**

pridemobility.com |    

# MORE ROOM TO BREATHE

Memory foam seat with 6" slider for adjustable support and premium, all-day comfort

Delta steering wheel with water-resistant LCD console

NFC-enabled key fob for enhanced security

Standard front basket and rear view mirrors

USB charger port

Enhanced LED headlight

Durable vinyl accents

Pneumatic tires mounted on split rims for off-road thrills

FDA Class II Medical Device<sup>1</sup>



## BIG PLANS REQUIRE THE BEST RIDE.

Durable and reliable, this mobility scooter offers supreme comfort thanks to full suspension and memory foam seating in several sizes. The innovative PX4 provides room to breathe – no matter your destination. Don't be afraid to think BIG!

### FEATURES

- NEW NFC key fob for enhanced security. Simply place or tap near the receiver to start
- 500 lb. weight capacity
- Impressive top speed of up to 9.3 mph
- Cruises up to 23.4 miles on a single battery charge
- Depth adjustable 18x18-20" memory foam seat with 6" slider with additional seat size options that include 20x18-20" and 22x20-22"
- Delta steering wheel with LCD console
- Standard mirrors and front basket to hold whatever you need for the road ahead
- Water-resistant user interface
- Full LED lighting package for maximum visibility
- Pneumatic tires mounted on split rims for excellent handling and traction on every terrain
- Supplemental hand brake for additional stopping power
- Choice of five unique vinyl accent color options with a 7-year outdoor durability rating

### AVAILABLE VINYL ACCENT COLORS



### ACCESSORIES

Only one rear accessory can be used at a time.

- 3" Single/Double Cane or Crutch Holder
- Wishbone Crutch Holder
- Cane Crutch Holder
- Forearm Crutch Holder
- Walker Holder
- Cell Phone Holder
- Oxygen Tank Holder
- Cup Holder
- Saddlebag
- Lap Belt
- Weather Cover
- Safety Flag
- Rear Basket

Cup holder, rear basket and saddlebag can be ordered as a bundle.


### MILITARY PATCHES



LCD digital display



Enhanced LED headlight



Full suspension