

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST COMFORTABLE RECLINER IN THE WORLD!

THE 4-MOTOR RECLINER



FDA Class II Medical Device¹

Shown in Canyon Walnut Fabric



New! Ergonomically designed toggle remote with two programmable memory positions



Three different heat pads located in shoulders, lumbar and seat



RADIANCE¹

Whether you cozy up with a beverage and your favorite book or kick back for a relaxing evening of television, this power recliner is right for you! Experience the warmth of the **Radiance Collection** from VivaLift!® Power Recliners. With three different heating pads in the shoulders, lumbar and seat, you can adjust your heat level accordingly and soothe your tired muscles. Keep your smart phone close at hand while you charge it with the wireless charging station. And if the power goes out, know you have peace of mind with the lithium battery backup.



Scan here for more information on **VivaLift!®**









PRIDE®
MOBILITY
LIVE YOUR BEST®

RADIANCE¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	RADIANCE ¹	RADIANCE ¹	RADIANCE ¹	RADIANCE ¹
MODEL	PLR-3955PW	PLR-3955S	PLR-3955M	PLR-3955LT
POSITION	True-Infinite	True-Infinite	True-Infinite	True-Infinite
SIZE	Petite Wide 5'4" and below	Small 5'4" and below	Medium 5'5"-6'1"	Large/Tall 6'2" and above
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	19.5"	18"	19.5"	21.5"
SEAT DEPTH	18"	18"	21"	22"
SEAT WIDTH	24"	20"	22.5"	23.5"
TOP OF BACK TO SEAT	23"	25"	28"	31"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	40"	41"	43"	49"
DISTANCE FROM WALL	14"	15"	15"	21.5"
OVERALL WIDTH	41"	37"	39"	40"
CHAIR WEIGHT	161 lbs.	156 lbs.	158 lbs.	178 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. True-infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

STANDARD FEATURES

7. NEW! Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions, plus turns your heat on/off and allows you to adjust heat levels for your comfort. Convenient USB charger port allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.


Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



ES QUICK SHIP = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

STANDARD FABRIC: CANYON



OCEAN

SILT

STEEL

WALNUT

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]