



## Transformer 2 Scooter

What are the dimensions unfolded? **37" high X 19 ½" wide X 39" long with tiller at highest**

What are the dimensions folded? **26" high X 19 ½" wide X 16 ½ long with tiller at lowest**

Lowest tiller height from floorboard? **19"**

Highest tiller height from floorboard? **31"**

What are body construction materials? **Aluminum Frame, Plastic Trim Pieces**

What is the height from top of the floorboard to ground? **6 ¾ "**

What is the seat thickness? **3"**

What is the seat made of? **Foam, Vinyl Covering**

What is the seat depth? (front to back) **15"**

What is the seat width? **17"**

What are the tires made of? **Solid Rubber**

What is the ground clearance (bottom of floorboard)? **3.25"**

Distance from front of seat to floorboard? **13"-16" seat cushion**

Distance from front of seat to tiller? **8"-14 ½ "**

How long does it take the battery to charge? **5-6 hours**

Can multiple remotes be programmed at the same time? **Yes, up to four remotes can be synced.**

Lowest Ground Clearance? **1" rear anti tipper wheels**